



NMR Operation at NYSBC

NYSBC Solid State NMR
Short Course

Schedule and outline, day 1

Bruker Acquisition and Processing Software
Navigating in TOPSPIN

Topspin Directory Structure

Data structure

Back-up/housekeeping protocols

Data server

Instrumentation

Magnet safety

Console components

Probes and probe safety

Initial setup

Samples and rotors: packing and spinning

Tuning and matching using Wobb

High power tuning and matching using oscilloscope

Variable temperature and coil air flow

Schedule and outline, day 2

Basic solid state NMR experiments

Setting up Magic Angle with KBr

Power handling

Direct excitation with CW decoupling

Composite pulse decoupling (CPD)

Cross Polarization (CP) with CPD

Shimming

More advanced solid state NMR experiments

Hahnecho

1D double CP experiment


^{13}C - ^{15}N correlation – 2D double CP experiment

^{13}C - ^{13}C correlation – RAD type polarization transfer

^{13}C - ^{13}C spin pair correlation – RFDR

^{13}C - ^{13}C - ^{15}N correlation – Doubcp-Radmix experiment

IT policies: spectrometer PC



Isolated from the outside world

Sftp out only

Raw and processed data deleted ~90 days

Parameter sets retained forever

IT policies: data backup

File not Found

Current user backed up to server every 10 min
Entire spectrometer backed up to a server once a week
Backed up to tape once a month
Intended for disaster recovery
Not intended for long term archiving

IT policies: data server




Net access to data

Sftp 'username'@sftp.nysbc.net

Replication of each spectrometer PC

Check on progress of data acquisition

Topspin directory structure



Twiki pages

<http://www.nysbc.net/twiki/view/Main>

<http://www.nysbc.net/twiki/register.html>

<http://www.nysbc.net/twiki/bin/view/Main/NmrfacGroup>

<http://www.nysbc.net/twiki/bin/view/Main/NmrSchedule>

<http://www.nysbc.net/twiki/bin/view/Main/CaL600>

Topspin



Starting topspin

Tab/window layout


Browser

Frequently used buttons

Manipulation of processed data

Miscellaneous features

Topspin interface



Topspin buttons

Buttons for vertical scaling (intensity manipulation)

A row of six buttons for intensity manipulation: *2, /2, *8, /8, a smooth icon (diamond with vertical lines), and a reset icon (square with a cross).

Increase the intensity by a factor of 2 [***2**]

Decrease the intensity by a factor of 2 [**/2**]

Increase the intensity by a factor of 8 [***8**]

Decrease the intensity by a factor of 8 [**/8**]

Increase/decrease the intensity smoothly

Reset the intensity [**.vr**]

Buttons for interactive manipulation

A row of seven buttons for interactive manipulation: phase (sine wave), integration (integral symbol), multi-display (three stacked bars), baseline (wavy line), peak picking (vertical line with a cross), calibration (triangle with a circle), and distance (arrow with a cross).

The functions of the individual buttons are:

Switch to phase correction mode

Switch to integration mode

Switch to multi-display mode

Switch to baseline correction mode

Switch to peak picking mode

Switch to calibration mode

Switch to distance measurement mode

Buttons for data handling:

A row of nine buttons for data handling: new (document with plus), open (folder), save (document with disk), email (document with envelope), print (printer), copy (document with arrows), paste (document with arrows), and 2d (2D icon) and 3d (3D icon).

The functions of the individual buttons are:

Create a new dataset [**Ctrl+n, new**]

Open a dataset [**Ctrl+o, open**]

Save the current dataset [**Ctrl+s, sav**]

Email the current dataset [**smail**]

Print the current dataset [**Ctrl+p, print**]

Copy the data path of the active data window to the clipboard [**copy**]

Paste the data path on the clipboard to the active data window [**paste**]

Switch to the last 2D dataset [**.2d**]

Switch to the last 3D dataset [**.3d**]

Buttons for display options

A row of four buttons for display options: units (h/p), y-axis (y-axis symbol), overview (spectrum icon), and grid (grid icon).

The functions of the individual buttons are:

Toggle between Hz and ppm axis units

Switch the y-axis display between abs/rel/off

Switch the overview spectrum on/off

Toggle grid between fixed/axis/off

NYSBC ultra high field magnets




800US#2_NE



800US#1_NW



900#2



800CONV



750WB



900#1

NYSBC high field magnets



500



600



700

Basic magnet safety

Magnetic objects

Keep ferrous/steel items outside 5G stray field (red line)

Gas tanks

Chairs

Electronics

Tools

Pens

Razor blades

Hair clips

Earrings

Medical implants

Keep medical implants outside 5G stray field (red line)




Magnet Quench



Evacuate the room immediately:
large volumes of N_2 and He may be released very quickly and
 O_2 concentration in the air will drop dramatically.
You will keel over and die.

Bruker Ultra Stabilized magnet



Bruker magnet PC

BMPC for ultra high field magnets (900s, 800s, 750wb)

Monitor magnet parameters:

P1, P2, P3, T1, T2, He flow

Control pumps for helium flow


Call out in event of problems




Spectrometer console



750wb spectrometer console



750wb spectrometer console




↑
Communication,
timing and
frequency control
units
(CCU, TCU, FCU)

↑
Receiver
(RCU)

↑
Signal generator
units (SGU 1-4)

↑
Routers

BLAX/BLAH amplifier



300W output
(HRMAS/solution)

Status LED lights

RF input,
1000W/300W select,
trigger

1000W output
(CPMAS)

Power supply
Voltage LED lights

EDASP

Main Frequency and Carrier offset

Channel

Signal Generator

Amplifier

Edit Spectrometer Parameter

frequency	logical channel	amplifier	preamplifier
BF1 188.642969 MHz	NUC1	<div> <div>X 1000 W</div> <div>X 300 W</div> <div>H/F 1000 W</div> <div>H 100 W</div> <div>X 1000 W</div> <div>X 300 W</div> </div>	<div> <div>1H LNA</div> <div>XBB31P 2HS</div> <div>2H</div> <div>HPHP 19F/1H</div> <div>HPHP XBB31P</div> </div>
SFO1 188.661833 MHz	F1		
OFS1 18864.3 Hz	13C		
FCU1/SQU1			
BF2 750.22 MHz	NUC2	<div> <div>X 1000 W</div> <div>X 300 W</div> <div>H/F 1000 W</div> <div>H 100 W</div> <div>X 1000 W</div> <div>X 300 W</div> </div>	<div> <div>1H LNA</div> <div>XBB31P 2HS</div> <div>2H</div> <div>HPHP 19F/1H</div> <div>HPHP XBB31P</div> </div>
SFO2 750.222 MHz	F2		
OFS2 2000.0 Hz	1H		
FCU2/SQU2			
BF3 750.22 MHz	NUC3	<div> <div>X 1000 W</div> <div>X 300 W</div> <div>H/F 1000 W</div> <div>H 100 W</div> <div>X 1000 W</div> <div>X 300 W</div> </div>	<div> <div>1H LNA</div> <div>XBB31P 2HS</div> <div>2H</div> <div>HPHP 19F/1H</div> <div>HPHP XBB31P</div> </div>
SFO3 750.22 MHz	F3		
OFS3 0.0 Hz	off		
FCU3/SQU3			
BF4 750.22 MHz	NUC4	<div> <div>X 1000 W</div> <div>X 300 W</div> <div>H/F 1000 W</div> <div>H 100 W</div> <div>X 1000 W</div> <div>X 300 W</div> </div>	<div> <div>1H LNA</div> <div>XBB31P 2HS</div> <div>2H</div> <div>HPHP 19F/1H</div> <div>HPHP XBB31P</div> </div>
SFO4 750.22 MHz	F4		
OFS4 0.0 Hz	off		
FCU4/SQU4			

• : cortab available

Save Switch F1/F2 Switch F1/F3 Default Cancel Param


1H high resolution (solution) preamp

0-31P high resolution (solution) heteronuclear preamp (2H bandstop)

19F/1H high power (solids) preamp

0-31P high power (solids) preamp (with sideways plugin bandpass filters)

750wb magnet space




AirJet VT line

Probe box

Trough with VT tubing
(next to the wall, behind the magnet)

Preamplifier (HPPR) tower

Preamplifier tower



Magic Angle Spinning

4mm ZrO rotors	2-15kHz
3.2mm ZrO rotors	4-23kHz
3.2mm Sapphire rotors	2-15kHz
1.3mm ZrO rotors	35-65kHz

Rotors should be well packed with a **homogeneous** sample.

Spacers should be **centered**.

Rotor caps should have **snug** fit and be in good **condition**.

A **black mark** should be drawn on the rotor bottom beveled sidewall. Make sure that the mark is fairly long. Some rotors have laser etched marks that are not bright enough – paint them over.

MAS controller is capable of starting spinning, maintaining spinning frequency within 2-5Hz, stopping spinning and ejecting a rotor.

Make sure that all air tubing quick connects are plugged in properly.
Make sure that MAS detector cable is plugged into the probe.

Rotor eject/insert on WB MAS probes

It is possible to eject and insert rotors in Bruker probes without removing a probe out of the magnet.

Before ejecting/inserting a rotor, **turn air flow** through the coil **off**.

Before ejecting a rotor, make sure that the transfer line is capped with a **safety cap**. However, don't close the safety cap all the way – leave a 1/8" opening to facilitate eject air flow.

To eject a rotor, press EJECT button on MAS controller. If a rotor doesn't eject, **gently and slightly** rotate ($<5^\circ$) and move up and down (2mm) the transfer line. Unplug and plug eject tubing quick connect into the transfer line.

To insert a rotor, make sure that there is no rotor already inside the probe. Turn the coil air flow off.

Drop a rotor into the transfer line with the cap pointing up.

Disconnect a 8mm gray tubing (shim stack purge line) from the quick connect joint on the shim stack brass ring and use this air flow to push rotor into the proper position inside the stator.

Start spinning before you turn air flow through the coil on.


If a rotor doesn't spin, eject and insert again.

MAS controller: front and back view




Comments:

When switching between SB and WB cable sets, don't pull air input cable.



Air input tubing WB air tubing set

MAS controller: front left




MAS controller: front right

Set spinning frequency

Start spinning

Adjust for rotor type and WB/SB



Solid state probes at NYSBC

Probe	Tuning range, MHz	RF fields, kHz	Comments
750MHz 4mm WB HXY	X:150-200 Y: 70-130	H/C/N: 100, 50, 40 DCP: 100/35/25	Good RF reliability. Unstable for VT DCP expts
750MHz 4mm WB HFX	H/F X: 40-220	H/X: 100, 50, 40	19F only 50kHz
750MHz 3.2mm WB HCN Efree and SB HCN	narrowband	H/C/N: 100, 60, 35 DCP: 100/35/25	The best probes to use for triple channel VT experiments
750MHz and 900 MHz 4mm HCND HRMAS	narrowband	H/C/N: 4us/6us/12us	
900MHz 4mm SB HCN Efree	narrowband	H/C/N: 100, 60, 35 DCP: 100/35/25	Good probe for triple channel VT expts
600MHz 3.2mm HCND and HCN DNP	narrowband	H/C/N: 100, 60, 35 DCP: 100/35/25	Testing/installation
900MHz 3.2mm SB HX	X:40-250	H/X: 80/50	Good for quadrupolar expts

750 WB HFX probe (1)




Comments:

Check that connections are tight and high power cables are screwed on all the way

Low Frequency (LF)
40MHz-220MHz channel

High frequency (HF)
H and/or F channel

750 WB HFX probe (2)



Comments:


Don't forget
frame cooling

900MHz spectrometer:
make sure that you
Use MAS cable,
not VT heater cable,
they have similar three
pin connectors

Frame cooling
quick connect

MAS cable
(three pins)

750 WB HFX probe (3)



VT heater

VT air line
connector


VT sensor

Comments:

Use the lower (TC2)
VT sensor

Use **utmost** caution with
glass ball joint
connector, do not crack
the glass dewar

750 WB HFX probe (4)



X channel tuning rod (TX)

Comments:

Before turning a tuning/matching rod make sure it's the correct one

Do not accidentally turn magic angle adjust (unless you are adjusting magic angle)

Initial setup

Samples and Rotors: packing and spinning

Tuning and Matching using wobb (low power tuning)


Tuning and Matching using oscilloscope (high power tuning)

Variable temperature and air flow


Bruker rotor/insert sets

Rotor type	Insert type	Active volume, μl	Active length, mm	B ₁ homogeneity	Shimming
Fully drilled	none	100	10	50%	Solids: 4/12Hz HRMAS: n/a
	tube	35	5.5	70%	Solids: 3/10Hz HRMAS: 1.5/15/25
1/2 drilled	large	12	1.7	90%	Solids: 1.5/4Hz HRMAS: 1.5/15/25
	small	38	???	>70%	Solids: 3/10Hz HRMAS: 1.5/15/25
2/3 drilled	small	60	7.5	>70%	The same
Fully drilled or 2/3 drilled	Custom Center	35-40	5	90%	The same
Fully drilled 3.2mm	none	40ul		90% for 18-20ul	Solids: 4/12Hz

Bruker rotor/insert sets: view



Rotor (un)packing tools



Packing rotors

Pack a rotor homogeneously.

Grind the sample uniformly

Center inserts

Don't leave empty spaces

Make sure that the cap:

Is in all the way inside the rotor

Fits snug inside the rotor

Is in good condition (check the fins with the magnifying glass)

Use only

ZrO (the same material as rotor) caps for high/low temperature experiments

Vestbule (brown plastic) caps for low temperature experiment.

tight fitting KelF (transparent plastic) caps for all other experiments

Power and voltage

Power (P) is proportional to voltage (V) and current (I) in the circuit:

$$P=VI$$

Accordingly to Ohm's law, where R is circuit resistance.

$$I=V/R$$

Substitution of the second equation into the first gives: $P=V^2/R$.

In an AC circuit, average voltage (V_{av}) is related to peak-to-peak voltage (V_{pp}) as: and complex impedance (Z) is used instead of resistance (R).

$$V_{av}=V_{pp}/2*\sqrt{2}$$

Therefore, in an AC circuit power is related to voltage as:

$$V_{pp}^2/8Z.$$

$$P = (V_{pp}/2*\sqrt{2})^2/Z =$$


Since impedance in all the NMR circuits is 50W,

$$P = V_{pp}^2/400.$$

Tuning and matching(1)

source and load

When a source sends power into a load,
some power is transmitted (dissipated in the load)
and some power is reflected back to the source




$$Z_{\text{source}} = Z_{\text{load}} \implies P_{\text{refl}} = 0$$

Tuning and matching(2)

source and load in NMR

In NMR, amplifier is a power source, probe is a load.
The B_1 field is generated by the power transmitted to the RF coil.
The objective is to minimize the reflected power.



Tuning and matching(3) NMR probe

Source impedance is always 50Ω

NMR probe impedance can be tuned and matched to 50Ω for a given frequency/range of frequencies

**We use tuning and matching adjustments:
capacitors, coils, variable length transmission lines**

Tuning and matching

If radiofrequency (rf) power is transmitted into a circuit, a part of it is reflected back

and a part of it is dissipated in the circuit.

To maximize the dissipation of the transmitted power in a probe, the probe impedance

at the rf frequency should be equal to the impedance of the outside circuit (50 Ohm).

It is achieved with tuning and matching capacitors.

The former adjusts the frequency of the circuit (tuning) and the latter optimizes the

power dissipation at this frequency (matching).

Rough tuning should be performed at low power (wobbling).

Precise tuning should be performed at high power (using oscilloscope).

Wobbing


**EDASP: make sure that EDASP pathway is correct.
Check the description in the manual.**

To wobb X or Y channel, use the corresponding preamp.


**To wobb ^1H , switch EDASP configuration to 100W ->
 ^1H LNA preamp and connect ^1H probe channel to LNA preamp.
After wobbing, change EDASP and hardware configuration back to
BLAH1000 -> directional coupler -> ^1H probe channel.**

Wobb


Wobb goes from low to high frequency: 15N -> 13C -> 1H
Commands: *wobb high*, *wobb f1*, *wobb f2*, *wobb ext50*



Wobb with preamp rack



High power tuning: setup



Power goes from transmitter to probe via a bi-directional coupler. The coupler splits off 10% of voltage (1% of power) passing through it. This voltage can be measured via oscilloscope. The ratio of reflected to forward voltage is the most precise indicator of how well a probe is tuned.

High power tuning: why

High power tuning is more precise than wobbling, particularly for **higher** frequencies.

If the isolation between channels is low ($<15\text{dB}$), the channels become coupled to each other and they have to be tuned simultaneously. It is possible only with the oscilloscope.

You should **always** monitor reflected voltage on the oscilloscope to make sure that no **arcing** is taking place.

If something is wrong with your hardware setup, parameter set, or pulse sequence, you can observe it on the oscilloscope.

High power tuning: how

Make sure that 20dB attenuators are plugged into the oscilloscope!
The oscilloscope input should **not** exceed 9V.

Use tuneHXY data set:
three channels pulsing simultaneously
d1=1s, aq=1ms, all pulses $\leq 1\text{ms}$, power levels 6dB


Measure forward peak-to-peak voltage

Measure reflected peak-to-peak voltage.

Minimize reflected voltage by adjusting tuning and matching knobs.
Overcouple to achieve better tuning/matching.

Aim for $<5\%$ reflected to forward voltage ratio.

Oscilloscope: front view



Oscilloscope: right side panel

Cursor menu
Cursor select
Cursor move

Quickmenu:
basic scope settings
(trigger, impedance)


Channel on

Trigger amplitude
level


Active channel off

Active channel
Amplitude scale

Active channel
Time scale




Oscilloscope: left side



Peak to Peak
voltage
between cursors

53
Trigger level

Oscilloscope: left side



Variable temperature on 750MHz spectrometer

WB and SB CPMAS Bruker probes are capable of cooling/heating a spinning rotor within +100/-100C temperature range.

A rotor is cooled/heated by a separate gas flow running through the center of the stator.

The gas flow is regulated by BVT3000 unit located in the spectrometer console.

A temperature sensor reads gas temperature inside the stator.

Use EDTE command to set gas flow and temperature from the console.

VT gas is delivered to the stator via an evacuated glass dewar.

Heating coil is located at the bottom of the dewar.

AirJet Cooler is used to cool the input gas (2cfm @ -85C input temperature).

Effective temperature inside a rotor depends on spinning frequency and gas flow.


NMR thermometer (lead nitrate or methanol) can be used to calibrate sensor and sample temperature.

Cooling setup at 750MHz spectrometer



Heater/controller

Cooler



Air input
from console

Manual
air flow
meter

Flow
regulator

STOP and READ this

If you have any questions/problems/concerns – STOP the experiment and call Boris at #207 or 917-526-1791.

Before you start acquisition (*go, gs, zg*), ALWAYS CHECK: acquisition time (*aq*), recycle delay (*d1*), *EDASP* configuration, and power levels.

Acquisition time default unit is SECONDS! Make sure that your acquisition time is ms.

Bruker power units are dB. The largest power setting is -6dB, the smallest is 120dB. Remember that 0dB is NOT zero power, it's 25-50% of maximum amplifier output.

Check power levels in the probe info page. Make sure that your powers and pulses are not longer than specified limits. If you have any doubts whatsoever, call Boris.

Always monitor reflected powers on oscilloscope. If you see arching, stop the experiment immediately.

Make sure that there is at least 1200 lph ambient air flow through VT line on HFX/HXY probes and 1600 lph on SBDVT probes.

STOP and READ some more

Always use the following configuration for CPMAS experiments:

SGU1 – BLAX1000 #1 – bidirectional coupler - high resolution broadband preamp -13C bandpass filter – probe for 13C

SGU1 – BLAX1000 #1 – bidirectional coupler - high power broadband preamp – filter – probe for any other nucleus except 1H and 15N

SGU2 – BLAH1000 – bidirectional coupler – 1H bandpass filter – probe for 1H (bypass the preamp if you are not observing 1H)

SGU3 – BLAX1000 #2 – bidirectional coupler - high power broadband preamp – 15N bandpass filter – probe for 15N

Always use the following configuration for HRMAS experiments:

SGU1 – BLAX300 #1 - high resolution broadband preamp - probe for 13C

SGU2 – BLAH150 -1H high resolution preamp – probe for 1H

SGU3 – BLAX1000 #2 - high power broadband preamp - probe for 15N

13C, 1H and 15N are cortab'd – the corresponding amplifier outputs are linear for this nuclei. However, the other nuclei are not cortab'd. The power output at the other nuclei frequencies will be slightly higher and non-linear (within 20%).

Use 1H high resolution preamp to wobb proton.

If you are running a VT experiment, make sure that shim stack/sample transfer line are purged properly (see VT section).

Magic angle



Spinning a sample at the magic angle $Q = 54.7$ ($1-3\cos^2Q = 0$) attenuates anisotropic interactions (CSA, dipolar, partially quadrupolar)

It is the only way to achieve high resolution spectra in solid samples

Setting up magic angle precisely and easily is very important

^{79}Br in KBr is convenient system to set a magic angle on

KBr is cheap and stable

^{79}Br has high natural abundance, and its gyromagnetic ratio is very close to that of ^{13}C

It is possible to set magic angle using KBr very quickly with 0.5° precision

For higher precision other systems should be used

Setting up magic angle on KBr

Load a KBr sample

Set up a zg experiment with ^{79}Br nucleus on channel 1
(or, better, pick it up from “bitin” directory)

Wobble:

if the tuning dip is anywhere close to the cursor, don't bother wobbling

Run a 1 scan zg experiment, ft and apk the spectrum,
set the O1 carrier exactly on resonance


Measure the first spinning side band/center band ratio

If it's less than 1:10 adjust magic angle
start acquiring in gs mode.

You should be able to see a train of rotationary echoes (next slide) on FID.


Very (!) slowly adjust magic angle to maximize rotational echoes.

Setting up magic angle on KBr: off magic angle



The carrier frequency should be set exactly on resonance to obtain FID as shown above. Clearly the first spinning side band/center band ratio is worse than 1:10 and the rotary echo pulse train is too short (it should last 3-4ms)
Slowly adjust magic angle rod to achieve this (next slide)

Setting up magic angle on KBr: on magic angle



Slowly adjust magic angle rod to achieve the FID shown above.
Now, the rotary echo pulse train extends for 3-4ms and
the first sideband/centerband ratio is 15%

Power handling of solid state probes

All solid state probes arc sooner or later..

No NMR spectroscopists are satisfied with probes' power handling capabilities..

Bruker WB HXY and HFX probes are spec'd for:

$^1\text{H}/^{13}\text{C}/^{15}\text{N}$: 120kHz@50ms, 50kHz@10ms, 50kHz@10ms in single channel mode.

We test the probes with a double cross polarization experiment:

$^1\text{H}/^{15}\text{N}$ 60kHz/50kHz@3ms followed by $^1\text{H}/^{13}\text{C}/^{15}\text{N}$ 100kHz/25kHz/35kHz@8ms

followed by ^1H 120kHz@1ms with 3s recycle delay

Do not exceed 80% of this total power input.

Always track reflected power on the oscilloscope.

Always test an adjusted/new pulse program on oscilloscope at 20dB lower power levels.

Make sure that your power levels are correct (remember 0dB is not 0% power!!!)

For longer decoupling times, decrease decoupling field proportionally

Power and voltage

Power (P) is proportional to voltage (V) and current (I) in the circuit:

Accordingly to Ohm's law, where R is circuit resistance.

Substitution of the second equation into the first gives:

In an AC circuit, average voltage (V_{av}) is related to peak-to-peak voltage (V_{pp}) as: and complex impedance (Z) is used instead of resistance (R).

Therefore, in an AC circuit power is related to voltage as

Since impedance in all the NMR circuits is 50W,

$$P = VI$$

$$I = \frac{V}{R}$$

$$P = \frac{V^2}{R}$$

$$V_{av} = \frac{V_{pp}}{2\sqrt{2}}$$

$$P = \left(\frac{V_{pp}}{2\sqrt{2}} \right)^2 \frac{1}{Z}$$

$$P = \frac{V_{pp}^2}{400}$$

Power, W	0.1	1	10	50	100	500	1000
Voltage, V_{pp}	6	20	60	140	200	450	630

Pulses, powers and dB units

Traditionally, attenuation of power and voltage in electrical circuits is measured in dB units.

We say that power is attenuated by N dB if it is decreased by a factor of $10^{N/10}$.

We say that voltage is attenuated by N dB if it is decreased by a factor of $10^{N/20}$.

The table below relates voltage and power attenuation for several dB values.

<i>dB</i>	<i>1</i>	<i>3</i>	<i>6</i>	<i>10</i>	<i>20</i>
<i>Power attenuation by</i>	0.26	2	4	10	100
<i>Voltage attenuation by</i>	0.12	0.4	2	3.16	10

Evidently, it is easiest to measure power in 3 or 10dB steps and voltage in 6 or 20dB steps.

Pulses, powers and fields in TOPSPIN

RF field is proportional to voltage and to square root of power.

$$B_1 \propto V, \quad B_1 \propto \sqrt{P}$$

Bruker software uses dB as a power unit.

-6dB - maximum amplifier output,

$$P_{\max}, \quad V_{\max}, \quad B_{1\max}$$

0dB – typically used amplifier output

$$\frac{P_{\max}}{4}, \quad \frac{V_{\max}}{2}, \quad \frac{B_{1\max}}{2}$$

120dB – zero power/voltage/rf field

Use command *pulse* to calculate powers/pulses/fields in Topspin:

Measure pulse/field at a certain power level

Switch to a zg data set (i.e. tune) and enter experimentally measured pulse/field values.

Use pulse command to calculate pulse length/field at a needed power level or vice versa

Composite pulse decoupling

There is never enough decoupling power in the solid state NMR of ^1H abundant samples.

Continuous wave (CW) decoupling is both less efficient and less broadband than composite pulse decoupling (CPD)

There is a number of CPD schemes in solid state NMR.

TPPM, XiX, and SPINAL are, probably, the most popular decoupling schemes.

The table below lists their advantages and drawbacks.

CPD	Advantages	Drawbacks
TPPM	The most efficient CPD for $<100\text{kHz}$ ^1H field	Pulse length and, sometimes, flip angle have to be optimized
SPINAL	The most efficient CPD for $>100\text{kHz}$ ^1H field	Less efficient than TPPM for low/medium ^1H fields
XiX	Fairly efficient for all decoupling fields. Doesn't have to be optimized	Less efficient than TPPM and SPINAL


Composite pulse decoupling (2)

TPPM




$P = 165\ (\theta), P' = 165\ (-\theta), \theta = 10^\circ\text{-}50^\circ$

SPINAL 64



$Q = 165(-10)\ 165(10)\ 165(15)\ 165(-15)X$
 $X165(20)\ 165(-20)\ 165(15)\ 165(-15)$


XiX



$P = 2.85 * \tau_r, \tau_r - \text{rotor period.}$

Lineshape broadening in SSNMR

Sample type	Interaction	What to do about it
Crystalline and ^1H abundant	^1H - ^1H homogeneous component in ^1H -X dipolar coupling	Decouple as hard as you can
Inhomogeneous and/or amorphous	Inhomogeneous line broadening	Nothing. Recrystallize it. Study another system.
Not ^1H abundant	Not ^1H - X dipolar coupling	Use medium strength ^1H decoupling field.
Paramagnetic/high unpaired electron density/etc	Very short T_2	Not much you can do about it.



Why CPMAS

CPMAS is the **most** often used 1D experiment in the world of biological NMR. Majority of ^{13}C and ^{15}N observe experiments start with a CP.

In a ^1H abundant spin system, $^1\text{H} - \text{X}$ CP increases sensitivity due to higher ^1H gyromagnetic ratio γ
shorter ^1H T_1 relaxation time.

CP is very easy to set up on a standard sample or on a system of interest.

Be aware

The spin system should be ^1H abundant and the $^1\text{H} - \text{X}$ distance should be fairly short

High molecular mobility may attenuate dipolar coupling to the point when polarization transfer rate becomes too slow.

Relaxation rate in spin locking frame $T_{1\rho}$ is generally fairly slow. However, presence of unpaired electrons or other relaxation pathways may substantially compromise CP efficiency.

Cross polarization


If in a IS spin system,
I magnetization is locked in the transverse plane with a field ω_I
and $\omega_I = \omega_S$ (Hartmann Hahn matching condition)

Then polarization is transferred from spin I to spin S at a rate proportional to
the effective IS dipolar coupling.
and magnetization enhancement equal to
 γ_I/γ_S

If a sample is spinning at ω_r spinning frequency, then Hartmann Hahn matching
condition changes to:
 $\omega_I = \omega_S + n \omega_r$, where $n = \pm 1, 2$


Also, in a cross polarization experiment, while S nucleus is observed,
relaxation time scale is determined by T_1 of I nucleus.

Hartmann Hahn matching – finger pattern



Hartmann Hahn matching condition in adamantane at SF = 5kHz.
Keeping ^{13}C field @ 50kHz and varying ^1H field in 0.25dB steps

Setting up CPMAS experiment



Pulse Sequence:	cp90.itin
SF (spinning freq)	10kHz
aq (acquisition time)	15ms
sw (sweep width)	350ppm
o1p (^{13}C carrier offset)	100ppm
o2p (^1H carrier offset, optimize)	5/+5ppm
d1 (recycle delay)	3-5s
pl1 (^{13}C CP field pl)	50kHz, 0/-3dB
pl2 (^1H CP field pl, optimize)	8/1dB
pl12 (^1H decoupling)	100k
(~0dB)	
p1/pl11 (pulse/power level to calibrate ^{13}C)	
P15 (spin locking time, optimize)	0.5-5ms

For spin lock, use 20-50% ramp for ^1H and maximum ^{13}C field.

Why glycine

While often enough it's possible to set up a CP experiment on the sample of interest, generally it's still better to use a model compound.

The model compound should be cheap, stable, reasonably well crystallized.

The model compound should have short enough (2-5 sec) T_1 and long enough ($>5\text{ms}$) $T_{1\rho}$.

The model compound should have
a CH_2 group (for decoupling optimization)
a $\text{C}(\text{O})$ group (for MA check/optimization)
an NH_3 group for N^{15} CP optimization.

n/a Glycine satisfies all these parameters.


Alternatives: NAV, other aminoacids.



Setting up CP: n/a glycine


C(O) peak:
lineshape indicates MA.
Expect 65-70Hz linewidth.

α -C peak. HH homonuclear
coupling makes CH₂ group
the most demanding for
decoupling. Expect 65-70Hz
linewidth



50ul centerpacked n/a glycine spectrum taken on 750MHz spectrometer.
CP with TPPM decoupling.
SF=10kHz, ns=4, ¹³C rf field 50kHz, ¹H rf field 120kHz. S/N = 170 (α -C).

Optimizing CPMAS: pop



#999, 998, 997, etc processed data

#899,898, 897, etc 2D data

Use pop macro to optimize the parameters:

Run a zg experiment and select a proper spectral window (right cursor button, dpl1)


Make sure that you are in processed experiment #1 and type pop command.

Enter a parameter you would like to optimize.

Decide whether or not you need to keep a 2D file

Start optimization.

Optimizing CPMAS: ^1H



Using popt to optimize ^1H parameters:

Set ^1H fields fairly low (<3 db) and measure p3 (1H 90).


It's more reliable to measure ^1H 180 – zero signal.

Once you know what ^1H field is, optimize pl2 (^1H spin locking field).

Pl2 optimal value may be different for C(O) and CH_2 .

Ramping ^1H field will shift optimal pl2 and make CP matching conditions less sensitive.


Optimizing CPMAS: more parameters



Measure ^{13}C field using p1/pl11
(optional ^{13}C 90 pulse)
Set the carrier on CH_2 peak (40ppm).
Set pl11 to approximately 50kHz (-3dB – 0dB)
Run popt with p1 varying from 1 to 10us
Zero signal (crossing) corresponds to 90 pulse.

Use popt to optimize TPPM decoupling.
Zoom in on CH_2 peak
Set pl12 to approximately 100kHz
(or lower for sensitive samples).
Run popt for p31 +/- 25% of ^1H 180 pulse.
Optionally: optimize o2 (^1H carrier offset)
in the range +/- 3000Hz.

CP polarization transfer dynamics



Optimal CP time depends on two parameters: polarization transfer rate and $T_{1\rho}$

Above, in a crystalline sample with reasonably slow relaxation rates (glycine), optimal polarization transfer time is 1ms for directly bonded ^{13}C (CH_2) and 6ms for 2 bond removed ^{13}C ($\text{C}(\text{O})$).

Optimizing CP on a real sample

Optimize CPMAS on glycine or NAV.

In a conductive sample, effective ^1H and ^{13}C fields may be significantly lower (by 10-30%) for the same power levels.

Nevertheless, **do not increase** power levels beyond recommended.

A good qualitative indication of sample conductivity is the ramping of reflected power. A conductive sample heats up during decoupling pulse and tuning/matching changes. It causes reflected power to change during the decoupling pulse – ramp.

If a sample is conductive, increase starting pulse widths somewhat and start optimization. Follow glycine optimization procedure.

If a sample is sensitive to heating (protein), start with **very low decoupling** power level (3-4dB down). You can always increase it later – better than baking a protein.

If a sample is mobile, keep spin locking time low – 0.5-1ms.

Use bigger ramp for large volume samples.

Optimizing CP on a real sample: more

If you can't see any signal with parameters fairly close to standard, quickly run pl2 optimization in 1dB steps.

If you still can't see any signal, try running direct observation with decoupling (hpdec).

If you still can't see any signal, put glycine back in and make sure that everything works properly.

If you can see hpdec signal, measure T_1 and T_2 .

If T_2 of the sample is very short (<5ms), the chances are that

T_{1r} is also too short for a successful polarization transfer.

Then, you will have to use hpdec.

Do you have protons in your sample?

Try ^1H direct observe just to make sure that you indeed have protons.

It's pointless to try $^1\text{H-X}$ CP in proton poor environment.

If you have a low density proton non-mobile sample with large $^1\text{H-X}$ distances, try long spin locking times – 5-10ms with large (50%) ramp.

Shimming

Typical NMR linewidth is 1-100Hz

Typical B_0 static magnetic field is 100-900MHz

Therefore B_0 field homogeneity should be 1ppb
(less than 0.1-5Hz)

Superconductive magnet coil has ~0.1ppt homogeneity in the center (sweet spot)

Extra homogeneity is achieved by:
superconductive/cryogenic shims
room temperature shims

Adjustable rf current runs through a 3D set of coils.
Resulting small magnetic fields are used to homogenize the main magnetic field

Superconductive shims are more powerful but they are not adjustable:
They bring linewidth down to 100Hz

A set of ~30 room temperature shims is used to bring linewidth down to 1-5Hz.

Shimming in SSNMR

In solid state NMR, shimming is performed by optimizing FID/FT signal lineshape of a rotating sample.

Typically, lineshape specs are
2-5Hz full width half height (FWHH) and 7-15Hz width at 10% height.

Generally, only few shims are needed to achieve lineshape specs.
Either ZY plane shims or ZX plane shims should be used.
Probes are usually designed so that the rotor is aligned in either ZY or ZX plane.

First order shims optimize FWHH: z, x²-y², z², xy(?) and
second order shims optimize a foot/shoulder of the signal lineshape: z², z³, z⁴

Shimming in SSNMR: a perfect sample

In the solid state NMR, shimming is performed by optimizing FID/FT signal lineshape of a rotating sample.

Typically, lineshape specs are 2-5Hz full width half height (FWHH).
FID is 200-500ms long, therefore decoupling has to be very low power (25kHz or less)

The sample should be highly ordered, have high signal (CP), short relaxation time and small HC effective coupling.

ADAMANTANE

Adamantane is a unique compound.


Molecular motion attenuates CH coupling dramatically.

25kHz TPPM decoupling is sufficient.

Adamantane is highly ordered and natural linewidth is less than 1Hz.

Adamantane has short (<3s) ^1H T_1 .

Shimming on adamantane: FID




This is an FID of a well shimmed small center packed rotor

The total FID length is 1s – that done just to show off, don't decouple more than 1/2s.

The beginning of FID looks fairly exponential – the foot will be small

Shimming on adamantane: FT, 60ul




This is more typical sample volume and shimming lineshape: 60ul

Full width half height:
z, x2-y2, z2 shims

The foot:
z2, z3, z4 shims
Use z4 to minimize asymmetric foot.
Use z3 to minimize symmetric linewidth at the bottom

Adamantane ^1H spectrum



In adamantane, ^1H homonuclear dipolar coupling is attenuated sufficiently to observe ^1H signal and check ^1H pulses/powers directly.

Imperfections



Nothing is perfect

**Pulse sequence commands take time
(typically few nanoseconds time scale)**


**Console electronic components need time to start
and stop pulsing**

**Hard rectangular pulses are not infinitely short.
Neither are they rectangular.
Rise and fall times are approximately 200ns.
Make sure that you take that into account when
measuring fields.**

**Long pulses are long and magnetization evolves
during them.**

**Amplifiers and probe need at least 4-5us to quiet
down after the last pulse before the acquisition
can start.**

B_1 homogeneity





In an infinite length solenoid coil B_1 is the same throughout the coil


In a limited length solenoid B_1 decreases towards the ends of the coil

Tapering a coil (tightening outside turns) may improve B_1 homogeneity somewhat

B₁ homogeneity




B₁ homogeneity: 2D ¹H field map



2D nutation experiment on a 60ul/5.5mm doped water sample.
The field map along rotor axis demonstrates 90% B₁ homogeneity for ~ 3.5mm length/35ul sample volume.
B₁ field drops off dramatically towards the ends of a rotor.

Amplifier linearity



SGU output is always linear

Amplifier output is non linear, particularly between 0 and -6dB

Cortab (correlation table) adjusts SGU output so that amplifier output is linear

Delays in pulse program and acquisition parameters

1u fq=cnst21:f2

10u pl12:f2 ;preselect pl2 drive power for F2

trigg

p3:f2 ph1 ;proton 90 pulse

0.3u

(p15 ph2):f1 (p15:spf0 pl2 ph10):f2 ;contact pulse

1u cpds2:f2 ;use cpdprg2=tppm15, SPINAL64 or XiX

d0

0.5u do:f2

(p1 pl11 ph3):f1 (1u do):f2

(mixing pl14 ph11):f2


(p1 ph5):f1 (1u cpds2):f2

Every frequency and power level switch takes time.

It is also safer to have some delays before and after composite pulses.

Always wait at least 4-5us before starting acquisition (dead time or d6).

Amplifier rise time




A typical Bruker or Varian (AMT) amplifier has approximately 200ns rise and fall times.

Therefore an experimentally measured 2.5us 90° pulse is, in fact, roughly 2.2-2.3us real pulse (corresponding to 110-114kHz instead of 100kHz field)

Use 360° pulse to measure field, not 180° – the relative error is smaller.

Ringdown



An NMR sample coil is not an ideal system.

A high power pulse generates a ringdown effect in the coil.

After the end of the pulse, residual currents keep running in the coil generating magnetic fields.


Ringdown time scale is proportional to probe efficiency (Q).

Ringdown time scale is inversely proportional to Larmor frequency.

Ringdown effect is usually
negligible above 300MHz: ringdown time < 5 μ s
very strong below 100MHz: ringdown time > 20 μ s

If signal lineshape is broad enough that 20 μ s of FID is important,
direct excitation signal will be distorted by ringdown effect.

Ringdown spectra





First several points of FID have abnormally high intensity – this is ringdown.

FT transform causes baseline roll.

Since the signal is fairly narrow, FID can be linearly predicted backwards or shifted to the left without any problem.

Convdt and left shift



If the ringdown time scale is much smaller than the FID time scale, it is possible to edit the beginning section of FID.

Use *convdta* command to convert Bruker digital FID into analog FID.

Set *datmod* = raw. Set *nsp* (number of set points) to the minimal number of points needed to eliminate ringdown

Use *ls* (left shift) command to shift FID to the left by *nsp* points. Process the resulting FID.

Linear Prediction (1)

Linear prediction (lp) is a method to construct FID points either forward or backward based on the existing FID section.

Lp forwards is used in 2D experiments to construct FID in the indirect dimension

Lp backward is used to counteract distorting effects of ringdown in the beginning of FID.

Linear prediction is tricky!

Linear Prediction (2)

For linear prediction backward:

Define ME_mod = lpr or lpc

Set TD_OFF to number of points in the beginning of FID you'd like to predict

Set TD_EFF to the number of points in FID you'd like to use for prediction

Set NCOFF to a “reasonable number”

For linear prediction forward for a indirect dimension 2D FID:



Define ME_mod = lfc or lfr

Set TD_EFF to the number of points you collected

Set LPBIN to the total number of points you would like FID to have

Set NCOFF to a “reasonable number”

Hahnecho (1)



If ringdown time is long enough on FID time scale, then neither left shift nor linear prediction would work reliably.


Hahnecho or solideocho refocusing scheme bypasses this problem.

A magnetization is allowed to evolve in transverse plane for either: rotor synchronized time period or as short time period as possible

Then, refocusing 180° (dipolar, hahnecho) or 90° (quadrupolar, solideocho) pulse is applied. After the same delay, the magnetization is more or less refocused. Since FID is collected after a substantial delay, ringdown is no longer an issue.

Homogeneous interactions (homonuclear dipolar coupling) are not refocused

Hahnecho (2)



hahnecho.itin

90 pulse p1


180 pulse p2

defocussing delay $n\tau_r$ d6


refocussing delay $n\tau_r - p2/2$ d7

Decoupling p31
pl12


T₂ measurement



$$I(t) = I_0 e^{-\frac{t}{T_2}}$$





Hahn-echo train:
Slow, needs good S/N,
precise





CPMG:
quick, works for low S/N,
has T1 error/contribution

T_1 measurement (inversion recovery)



$$I(t) = I_0 \left(1 - e^{-\frac{t}{T_1}} \right)$$

Measuring T_1 and T_2



Hahn echo can be used as a pseudo 2D experiment to measure T_2 . The signal intensity can be fit as a function of n (remember: rotor synchronized)

Inversion recovery experiment is used to measure T_1 . D array consists of exponentially increasing delay values D_i i.e. (1ms, 4ms, 16ms, 64ms, ...). Signal intensity can be fit as a function of delay

Double Cross Polarization

Double cross polarization is one of the most widely used heteronuclear correlation experiments in solid state NMR.

DCP is robust, easy to set up and effective way to achieve heteronuclear polarization transfer.

Typical DCP experiment used in an HCN spin system consists of:

$^1\text{H}/^{15}\text{N}$ cross polarization

followed by evolution t_1 time period to evolve ^{15}N chemical shift term,


followed by $^{15}\text{N}/^{13}\text{C}$ cross polarization,

followed by ^{13}C acquisition.


1D DCP experiment on a model compound can be used to optimize all power/pulse parameters.

U-gly (or another amino acid) is a good compound for DCP optimization.

2D DCP pulse sequence



Setting up DCP



Pulse Sequence:	cp90.itin
SF (spinning freq)	10kHz
aq (acquisition time)	15ms
sw (sweep width)	350ppm
o1p (^{13}C carrier offset)	100ppm
o2p (^1H carrier offset, optimize)	-5/+5ppm
d1 (recycle delay)	3-5s
pl1 (^{13}C CP field pl)	50kHz, 0/-3dB
pl2 (^1H CP field pl, optimize)	8/1dB
pl12 (^1H decoupling)	100kHz (~0dB)
p1/pl11 (pulse/pl to calibrate ^{13}C)	107

DCP: optimizing parameters

DCP is a very demanding pulse sequence in terms of power handling.

Start with a model compound (u-gly)

Optimize $^1\text{H}/^{15}\text{N}$ cross polarization, and measure $^{15}\text{N}/^{13}\text{C}$ fields using CP90 experiment.

Use ^1H decoupling field during $^{15}\text{N}/^{13}\text{C}$ CP such that $\omega_{\text{H}} > 2.5\omega_{\text{C}}$
Otherwise, polarization leaks from ^{13}C spin bath into ^1H spin bath.


Make sure that $\omega_{\text{C}} = n\omega_{\text{r}}$, $n=1, 2, 3(?)$
To achieve higher efficiency of $^{15}\text{N}/^{13}\text{C}$ polarization transfer

Within the limits described above, choose proper fields for $^1\text{H}/^{15}\text{N}/^{13}\text{C}$ section
(for instance 80kHz/35kHz/25kHz at 10kHz spinning frequency)

Use 10-20% linear or tangential ramp for DCP.

For a set ^{13}C DCP power level, optimize ^{15}N pl starting with calculated values.

Ubiquitin 2D DCP spectrum



Ubiquitin 50ul, 4mm rotor, HXY probe, 750MHz spectrometer,
75kHz ^1H decoupling, 35kHz ^{15}N field, 25kHz ^{13}C field,
ns = 16, t=-30C

^{13}C - ^{13}C correlation – RAD/DARR/PDSD


There are numerous ^{13}C - ^{13}C correlation methods in solids.
We are not going to discuss them here.

RAD/DARR/PDSD set of experiments

easy to implement
not hardware demanding
Work well in a u-labeled system
Recover all interactions within 4+ Å

non-specific
can't measure distance


RAD/DARR pulse sequence



Radmix.lorieau

13C 90 pulse	p1
13C 90 pl	pl14
1H 90 pulse	p3
1H 90 pl	pl3
Spin locking time	p15
Spin locking 1H pl	pl2
Spin locking 13C pl	pl1
Low power CW	
spin diffusion field	pl11
Low power CW	
spin diffusion delay	d9
TPPM decoupling (=SF)	pl12

Pf1 DARR spectrum



Pf1, 4mm rotor, HXY probe, 750MHz spectrometer,
80kHz 1H decoupling, 50kHz ^{13}C field

3D DCP-DARR pulse sequence

